



ThermoFisher
S C I E N T I F I C

Attune® NxT™ Acoustic Focusing Cytometer

Two Day Basic Training

Table of Contents – hyperlinked

[General Information/Support](#)

Section 1: Instrument Overview

[Fluidics](#)

[Optics](#)

[Electronics](#)

Section 2: Workflow

[Instrument start up](#)

[Performance test](#)

[Experiment setup & optimization](#)

[Compensation](#)

Data Acquisition

Data Analysis

[Instrument shutdown](#)

Section 3: [Maintenance](#)

Section 4: [General Software Overview](#)

Agenda – Day 1

9:00 am – 12:00 pm

Welcome and Overview of Training Agenda

Introduction to Flow Cytometry

Optional - Activity: Understanding Scatter

Attune® NxT™ Systems – Fluidics and Optics

Optional - Activity: Know Your Optics System

Attune® NxT™ Systems - Electronics

Software Overview & Experimental Workflow

12:00 pm - 1:00 pm

LUNCH

1:00 pm – 4:00 pm

Startup and Performance Tracking

Maintenance Overview

Voltage Walk

Agenda – Day 2

9:00am – 12:00 pm

Review Day 1

Attune® NxT™ Maintenance, Data & Account Management

Optimizing Instrument Settings

Instrument Startup and Performance testing

12:00 pm – 1:00 pm

Lunch

1:00 pm – 4:00 pm

Customer Samples and Multicolor Compensation (if applicable)

Autosampler (if applicable)

Review of Training – Additional time for practice

Instrument Shutdown


Hands-on Training Checklist

Attune Training Checklist			
<input type="checkbox"/>	Startup	<input type="checkbox"/>	How to get service
<input type="checkbox"/>	Shutdown	<input type="checkbox"/>	Advance settings - ASF
<input type="checkbox"/>	Performance test/baseline/LevyJennings	<input type="checkbox"/>	How to get follow up support
<u>Fluidics Functions/Maintenance</u>		<u>Software</u>	
<input type="checkbox"/>	Rinse	<input type="checkbox"/>	Create experiments/tubes
<input type="checkbox"/>	SIP Sanitize tube	<input type="checkbox"/>	Create/use/manage templates
<input type="checkbox"/>	SIP Sanitize plate	<input type="checkbox"/>	Parameters
<input type="checkbox"/>	Deep clean	<input type="checkbox"/>	Voltages/voltage walk
<input type="checkbox"/>	Sample Recovery	<input type="checkbox"/>	Threshold
<input type="checkbox"/>	Unclog	<input type="checkbox"/>	Workspace/Edit gates
<input type="checkbox"/>	Debubble	<input type="checkbox"/>	Create experiments/plates
<input type="checkbox"/>	System decontamination	<input type="checkbox"/>	Collection panel/run protocol
<input type="checkbox"/>	Correct way to remove/refill fluid bottles	<input type="checkbox"/>	Customize panel
<input type="checkbox"/>	AAS plate calibration	<input type="checkbox"/>	Compensation
		<input type="checkbox"/>	Options Menu
		<input type="checkbox"/>	Heat Map and HM set up
<u>Data Analysis</u>		<u>Data Management</u>	
<input type="checkbox"/>	Plots/gates	<input type="checkbox"/>	Export (exp/statistics/fcs files)
<input type="checkbox"/>	Statistics	<input type="checkbox"/>	Import (exp, fcs files)
<input type="checkbox"/>	Results tables	<input type="checkbox"/>	Data Base Management (Admin)
<input type="checkbox"/>	Overlaps		

<http://www.thermofisher.com/us/en/home/life-science/cell-analysis/flow-cytometry/flow-cytometry-learning-center.html>

Home > Life Sciences > Cell Analysis > Flow Cytometry > Flow Cytometry Learning Center

Flow Cytometry Learning Center



Learn about flow cytometry methods and technologies

The purpose of this learning center is to connect scientists (whether new or experienced) to our many resources for learning about flow cytometry applications, techniques and basic principles by providing a few key points of entry into the vast content.



Flow Cytometry Guided Learning

Find courses for basic applications and techniques used in flow cytometry

Molecular Probes School of Fluorescence

- [Flow Cytometry Basics](#)
- [Fluorescence Basics](#)

eLearning Courses

- [T Cell Stimulation and Proliferation](#)
- [Antibody Validation*](#)



Flow Cytometry Subtopics

Find webinars, application notes, white papers, videos and more for key flow cytometry applications and techniques

- [Areas of Biology](#)
- [Cell Isolation, Expansion and Differentiation](#)
- [Controls, Compensation and Calibration](#)
- [Flow Cytometry Assays](#)
- [Flow Cytometry Instrumentation](#)
- [Immunophenotyping](#)
- [Panel Design and Multicolor Flow Cytometry](#)

Web Technical Resources Continued

www.thermofisher.com/flowcytometry

ThermoFisher
SCIENTIFIC

Search All

Search



Contact Us

Sign In ▾

Quick Order



Flow cytometry

We are committed to accelerating your science by providing a comprehensive suite of solutions for the analysis of cells and their function by delivering you our flagship flow cytometry products designed to deliver high-performance results and save you time.



Flow cytometry antibodies, assays and reagents

Our diverse collection of antibodies, assays, beads and buffers support leading research areas.

- [Antibodies](#)
- [Assays and reagents](#)
- [Beads for instrument controls and standards](#)
- [Sample preparation buffers and reagents](#)
- [Custom antibody services](#)



Flow cytometry instrumentation

Our instruments are designed to make flow cytometry available to both new and experienced researchers.

- [Attune NxT Flow Cytometer](#)
- [Consumables and accessories for the Attune NxT Flow Cytometer](#)
- [Robotic automation for Flow Cytometry](#)



Flow cytometry education and support

Find tools, protocols and more detailed information or gain a basic understanding of techniques to help you [plan and execute your experiments](#).

- [Flow Cytometry Panel Builder^{beta}](#)
- [Flow Cytometry Protocols](#)
- [Invitrogen eBioscience Resources](#)
- [Flow Cytometry Learning Center](#)
- [Flow Cytometry Support Center](#)

Flow Cytometry Panel Builder^{beta}



This online tool guides you through flow cytometry panel design, providing a simplified, customizable experience to fit your flow cytometry panel design needs.

Molecular Probes School of Fluorescence



Fluorescence education for scientists. Learn the basics of fluorescence, imaging and flow cytometry through our Invitrogen Molecular Probes School of Fluorescence.

New →

Need additional training?



Cat #	Customer site
A25623	2 day NxT training
4484605	1 day NxT extended training
TRN00045	Hourly charge consulting /training

Detailed descriptions: <https://learn.thermofisher.com/flowcytometry>

For more information, questions or quotes - contact your FAS

www.thermofisher.com/aspire

Enroll now

Get your free product, rewards, benefits, and more

Rediscover your scientific aspirations


aspire
member program

Enroll now

Get your free product, rewards, benefits, and more

Join the free member program that supports your lifelong love of science.


Enroll today and receive a free, full-size product, 500 points towards rewards and discounts, and access to science career benefits and more...



Free full-size product

Select one free product from our catalog of nearly 20,000 kits, reagents, antibodies, and more.


Learn more >



Rewards and discounts

Easily earn points to redeem for rewards—apparel, electronics, lab items, and more.

Learn more >



Science career accelerators

Access valuable educational resources and benefits to help build your career in science.

Learn more >

Flow	4484604	FLO-103: Flow Cytometry FAS On-Site Training (1 Day)	Flow Cytometry On-Site Training	90,000	1
------	-------------------------	------------------------------------------------------	---------------------------------	--------	---

9

ThermoFisher
SCIENTIFIC

Obtaining Support

1. Check the troubleshooting section in the users guide.

If additional help is required, don't hesitate to contact: **1-800-955-6288**

(9:00am-8:00pm Eastern, 6:00am-5:00pm Pacific)

You will need the S/N of your instrument

Choose **Option 3** for Instrument Service and HW support

then **option 1** for a) instrument service or

b) to check status of instrument repair

Option 4, then **option 6**

Fluorescent Labeling/ Detection and

Flow Cytometry Reagents and Flow Applications

3. Or Email:

flowsupport@thermofisher.com

or

Instrumenthardwaresupport@thermofisher.com for TAC/Instrument services

Web submission: <https://www.thermofisher.com/us/en/home/support/instrument-repair-request.html>

Introduction to Flow Cytometry

What is Flow Cytometry?

CYTOMETRY is the measurement of physical or chemical characteristics of cells or particles

FLOW CYTOMETRY measurements are made as individual cells or particles in flowing stream pass through a flow cytometry instrument

- Performed on single cell suspensions
- Provides discrete measurements from each cell in the sample
- Provides a distribution of the measured characteristics in the sample

What makes a Flow Cytometer?

Abbreviated : FCM

Flow Cytometer is made up of 3 subsystems:

- **Fluidics**

To introduce and focus the cells for interrogation

- **Optics**

To generate and collect the light signals

- **Electronics**

To convert the optical signals to proportional electronic signals for computer analysis





Fluidics

Optics

Electronics

The purpose of a fluidics system is to transport particles in a fluid stream to the laser beam for interrogation


For optimal illumination, the stream transporting the particles should be in the center of the laser beam

Only one particle should move through the laser beam at a time

Fluidics system needs to be free of air bubbles & debris

Traditional Hydrodynamic Focusing

Particle positioning in laser is important




Acoustic Focusing

High sample input flow rates allow for more sample flexibility


12.5 $\mu\text{L}/\text{min}$

1,000 $\mu\text{L}/\text{min}$

Narrow particle focus = Narrow distribution




Narrow particle focus = Narrow distribution




Prior to wrapping in sheath


End-on view of capillary




Acoustic Focusing




Acoustic focusing off



Acoustic focusing on



Acoustic focusing on
(dilute sample)





Dilution
is the
solution!

Acoustic Focusing Capillary

- Acoustic Waves – similar to ultrasound used to visualize a fetus *in utero*
- 2.5 MHz

~20cm



Concentration and Flow Rates

The event rate will approach maximums stated in the column header when samples of stated concentrations are run at the flow rates below.


When acquiring large event files (i.e files with $> 10^6$ events), plot parameters should not be changed while recording.

Sample flow rate	<u>Fastest</u> (35,000 ev/sec) <i>maximum sample concentration</i>	<u>Accurate counts</u> (8,000 ev/sec) maximum sample concentration	Cell size and flow rate recommendation
1000 $\mu\text{L}/\text{minute}$	2.1×10^6 cells/mL	0.48×10^6 cells/mL	<ul style="list-style-type: none"> - Particles $> 4 \mu\text{m}$ - Predominantly acoustic focusing
500 $\mu\text{L}/\text{minute}$	4.2×10^6 cells/mL	0.96×10^6 cells/mL	<ul style="list-style-type: none"> - Particles $> 2 \mu\text{m}$ - Predominantly acoustic focusing
200 $\mu\text{L}/\text{minute}$	6.7×10^6 cells/mL	1.5×10^6 cells/mL	
100 $\mu\text{L}/\text{minute}$	1.3×10^7 cells/mL	3×10^6 cells/mL	
25 $\mu\text{L}/\text{minute}$	5.4×10^7 cells/mL	1.2×10^7 cells/mL	<ul style="list-style-type: none"> - Small particles $< 2 \mu\text{m}$ - Predominantly hydrodynamic focusing - Smallest sample core - Best resolution from background for dimly positives assays
12.5 $\mu\text{L}/\text{minute}$	1.0×10^8 cells/mL	2.4×10^7 cells/mL	

Let your biology and data quality be your guide. If good data is obtained while running at 2-8,000 ev/sec, adjust the sample concentration and flow rate to maintain that.

Attune® NxT Fluidics System

Focusing Fluid Reservoir




Valve

1ml Sample Syringe

Fluid Lines

Focusing Fluid Filters

Sensor Connections




Waste

Focusing Fluid

Wash Fluid

Shutdown Fluid

Attune NxT Fluidics System



Attune NxT Fluidics Solutions

Attune® Focusing Fluid: a buffered, azide-free solution which transports focused particles to the flow cell for laser interrogation. It prevents sample from coming into contact with the walls of the flow cell. It contains an anti-microbial agent, a preservative and a detergent designed to minimize bubble formation. Fluid must be RT before use.

Attune® Wash Solution: a solution for removing cellular debris and dyes from the fluidic system of the instrument.

Attune® Shutdown Solution: a solution which minimizes bubble formation and crystal deposit in the fluidics system when the instrument is shutdown.

Attune® Debubble solution: a solution formulated to remove bubbles from the fluidics system.

NEW! Attune® Flow Cell Cleaning solution: a solution which when diluted, will minimize contamination buildup that may occur in system lines or the flow cell and system lines.

10% Bleach: Fresh solution used to decontaminate the fluidics lines.

Deionized water: Used for diluting bleach. High quality, filtered and sterile.

About Bleach

10% bleach is defined as a 1 in 10 dilution (1 part bleach to 9 parts water) of 5.25% sodium hypochlorite in water. This gives a final concentration of 0.5% sodium hypochlorite equivalent to 5000 ppm of available chlorine.

Recently more concentrated formulations (Ultra and Concentrate) have become available:

Ultra is 6.15% Sodium Hypochlorite and should be diluted 1 part bleach to 11 parts water.

Concentrate is 8.25% Sodium Hypochlorite and should be diluted 1 part bleach to 15 parts water


Bleach Solution	Dilution	Chlorine (ppm)
5.25%	None	52,500
	1:10	5,250
Ultra 6.15%	None	61,500
	1:12	5,125
Concentrate 8.25%	None	82,500
	1:16	5,150

http://www.cdc.gov/hicpac/Disinfection_Sterilization/19_00glossary.html

Recommendation:

- Prepare fresh bleach
- Use laboratory-grade bleach
- DO NOT USE bleach with additives (such as perfumes or soap)

Fluidics Functions (on the Instrument tab)



- **Stop** - click to end any running routine.
- **Recover Sample** - returns unused sample volume back to the well or the tube.
- **Rinse** - flushes system between samples to minimize carryover. Rinse runs automatically every time the SIP is lowered, but it can also be user-initiated.
- **Sanitize Attune SIP** - sanitizes the SIP and sample lines between sticky/dirty samples or experiments; requires 10% bleach solution.
- **Deep Clean** - thoroughly washes the system sample lines and flow cell between sticky/dirty samples or experiments; requires 10% bleach solution (can also use debubble solution)
- **Startup** - primes the instrument fluidics with Attune® Focusing Fluid.
- **Shutdown** - automatically cleans, sanitizes and shuts down the instrument.
- **Debubble** - clears bubbles from the fluidics lines of the cytometer; Attune® Debubble solution required.
- **Unclog** - back flush operation to remove clogs from the sample line.
- **Decontaminate System** - semi-automated decontamination of the Cytometer and the Auto Sampler fluidics.
- **New! Export logs for Service** – one button export of log files for field service assistance

New – Attune NxT Flow Cell Cleaning Solution

What is Attune NxT Flow Cell Cleaning Solution?

Alkaline liquid concentrate

Removes contaminants without damaging the quartz flow cell

- Sticky Cells
- Dyes like PI


How often should it be used?

Daily - high volume users – running >6 hrs or 8 plates/day

Weekly – lower volume users


Directions

- 1) Combine 1 mL of Cleaning solution with 2 mL ultra pure water
- 2) Run Sanitize SIP (substituting the cleaning solution for bleach)
- 3) Prepare fresh 10% bleach
- 4) Run Shutdown




Cleaning Functions

Sanitize SIP


-  Quick wash/sanitize of sample line.
- Duration: 1 min. cycle time.
- Requires 3 mL 10% Bleach.
- Run between users especially after sticky samples, DNA stains or beads.

Deep Clean


- Sanitize system with bleach and wash solutions for selectable period of time.
-  Three levels: Quick (10 cycles/ 25 min), Standard (20 cycles/50 min, Thorough 30 cycles/75 min.

Shutdown

System clean and flush with bleach, wash, and shutdown.

- 

	Quick	Standard	Thorough
	10 cycles /30 minutes	20 cycles/60 minutes	30 cycles/75 minutes
- Few samples
- Immunophenotyping, apoptosis
- Samples with sticky dyes (PI), NLNW



Instrument placed in stand-by (dream state) upon completion


Located on the instrument tab

Instrument Cleaning Guide

Between samples	<ul style="list-style-type: none"> • Rinse – automatically initiated when SIP is lowered (for tubes), or set in <i>run protocol</i> for plates • Sanitize SIP between sticky samples or cell counts
Between users / experiments USE: 1) if there is ≥ 30 min between users. 2) If there is < 30 min between users.	1) Unclog then Quick Deep Clean - 30 minute cleaning routine (click on the arrow below the Deep Clean icon to select <u>Quick</u>) or 2) Unclog then 2X Sanitize SIP / Sanitize Autosampler SIP (plate experiments) – 1 st time with 3 mL 10% Bleach 2 nd time with 3 ml Wash or De-bubble solutions
End of day (3 steps)	<ul style="list-style-type: none"> • <i>Unclog</i> • <i>**SIP Sanitize with 1:3 dilution of Attune Flow Cell Cleaning solution</i> • <i>Thorough Shutdown</i> (click on the arrow below the Shutdown icon to select <i>Thorough</i>)

Note: Always wipe the outside of the SIP after doing a SIP Sanitize

Sample Recovery - Tubes



When is **sample recovery** available?

Anytime sample remains in the sample loop
Stop option has been reached
Operator clicks stop


NOTE: Must select **recover sample** before lowering the tube lifter

Sample Recovery - Plates


Sample recovery is disabled when moving between wells except when acquisition is stopped.

When Stop is pressed and the next well has been pre-loaded, sample recovery will recover the preloaded sample back into the plate.

1. Press the “Recover Sample” button




2. If sample has been preloaded, follow instruction to recover sample into plate
3. The remaining sample will be dispensed. This includes the dead volume and sample that has not been acquired



Status Indicator Lights



LED Color	Function/step
Blue (fade)	Warm up
Blue (solid)	Warm up complete
Blue (flashing)	Startup and instrument functions (except rinse)
Green (solid)	Startup complete or instrument idle
Green (flash)	Data/sample acquisition
Multicolor (fade)	Shutdown complete/sleep mode
Amber (blink)	Instrument Errors




Fade: Warming Up
 Flashing: Startup and instrument functions
 Solid: Warm up complete



Flashing: Instrument Error



Solid: Start up complete or Idle
 Flashing: Data/Sample Acquisition



Fade: Shutdown Complete/Sleep Mode

Fluidics


Optics

Electronics



Optical System Anatomy

- Up to 4 lasers
- Up to 14 fluorescence detectors
- 2 scatter detectors
- Different PMTs for different Wavelengths



Attune NxT Optical Components

- From 1 to 4 Lasers

Violet	405 nm	50 mW
Blue	488 nm	50 mW
Green	532 nm	100 mW
Yellow	561 nm	50 mW
Red	637 nm	100 mW


- Option of Green OR Yellow laser

- Filters are user changeable


- 2 Attune® NxT Accessory Filter Configurations

Violet Side Scatter kit

Fluorescent Protein Optimization kit




Flat Top Laser Profile



- Flattened gaussian profile
- Stable laser alignment
- Minimized down time
- Corrects for fluidic instabilities




Visible Light



What happens to light when it hits a cell?

Laser Light Scatter

- When laser light interacts with a cell, light is scattered in all directions
- We look at Forward Light Scatter and Side Light Scatter




Scatter


Laser Light Scatter

Forward Scattered light (FSC) is relatively proportional to cell-surface area or size

Side-scattered light (SSC) is relatively proportional to cell granularity/internal complexity of the cell. SSC is usually collected at 90 degrees to the laser beam




Flow Cytometry – FSC (Size)



<https://www.thermofisher.com/us/en/home/support/tutorials.html>


Flow Cytometry – FSC (Size), SSC (Complexity)




Laser Light Scatter

Measurements of FS and SS

- allow for differentiation of cell types in a heterogeneous cell population
- look for changes in cell health




Ammonium chloride lysed whole blood



Aged culture of Jurkat T cells: green are live cells & red are dead cells


Fluorescence

1.




Absorption

2.




Excitation




Light measured at 90°



3.



Emission



Flow Cytometry - Fluorescence




<https://www.thermofisher.com/us/en/home/support/tutorials.html>

Fluorescent Light – Common Definitions


- **Absorption (Excitation) spectrum:** The wavelength range over which a fluorescent compound can be excited
- **Emission spectrum:** The range of emitted wavelengths of a fluorescent compound, it is a longer wavelength than the absorption wavelength

Summary




- **Auto-Fluorescence:** fluorescence that originates from endogenous sample constituents which are excited by lasers and detected on PMT usually between 400-600nm.
- **Non-specific fluorescence:** Unbound or nonspecifically bound probes – this may increase 'apparent' autofluorescence.

Filtered Light Emission Path – Generalized Configuration




Optical filters and mirrors to direct and split light




Filters


Bandpass Filter




Longpass Filter




Shortpass Filter




Dichroic Mirror




Filtering Emitted Light



Blue - Standard Configuration




Violet Blue Yellow Red - Standard Configuration




Blue and Violet Side Scatter (with NLNW filter kit installed)


Filter Configuration



No-Wash, No-Lyse Detection of Leukocytes




Violet SSC Configuration





Violet Side Scatter Kit

415 dichroic LP
405/10 BP filter
Blank 25 mm holder



Fluorescent Protein Configuration

mCherry 615/25 BP filter
 GFP 510/10 BP filter
 YFP 540/30 BP filter
 525 dichroic LP



Fluidics
Optics
Electronics



Functions of Electronics:


- Converts detected light signals into proportional electronic signals (voltage pulses)
- Electronic signals are processed by the on-board processor
- Converts electronic signals from the detectors into digital data used for analysis
- Interface with the computer for data transfer

Flow Cytometry Detectors




Photomultiplier tube (PMT):




- Often referred as the “detector”
- PMT convert photons into electrons and amplify them to create a voltage pulse.




Sample Presentation: Voltage Pulse



Sample Presentation: Voltage Pulse



Pulse to Parameters






<input checked="" type="checkbox"/>	Time							
<input type="checkbox"/>	Event count							
Enabled			A	H	W			
<input type="checkbox"/>	Target		Label			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input checked="" type="checkbox"/>	FSC	<input type="text"/>	<input type="text"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input checked="" type="checkbox"/>	SSC	<input type="text"/>	<input type="text"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input checked="" type="checkbox"/>	BL1	<input type="text"/>	<input type="text"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input checked="" type="checkbox"/>	BL2	<input type="text"/>	<input type="text"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input checked="" type="checkbox"/>	BL3	<input type="text"/>	<input type="text"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input checked="" type="checkbox"/>	RL1	<input type="text"/>	<input type="text"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	RL2	<input type="text"/>	<input type="text"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	RL3	<input type="text"/>	<input type="text"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Attune default: All parameters (A-H-W) on all channels

Doublet Discrimination – 3 Ways to Display (blood)

In most cases, data analysis should include gating on single cells.



Pulse height = pulse area




Attune® NxT AutoSampler

Attune® Autosampler

- Compatible with many different standard plate formats, including 96-well, flat, round and V-bottom.
- Round Bottom plates are recommended for optimal mixing
- Intelligent probe design minimizes clogging and carryover (<0.5%) and prevents damage to the instrument
- Performs automated cleaning between wells (from 1 to 10 rinses) and when the instrument is shutting down
- Minimal variation regardless of sampling method (tube vs. plate) and collection rates
- Easy to plug and unplug on the Attune® NxT™ Cytometer



Attune® Autosampler Mixing procedure



The user sets:


- The plate type
- The total sample volume
- The number of mixes
(Max. 3 mixes, recommended at this time)

The system defines:


- The liquid level in well
- The probe position
- The mixing method

Mixing sample by aspiration instead of shakings ensures homogeneity of the sample and maintains cell viability

Exterior Components




Power Switch RS-232



Interior Components



Sample Probe



Leak Sensor



1ml Syringe

Waste

Focusing
Fluid

Instrument Start-up

Performance test

Create experiment

Experiment settings optimization

Compensation

Data acquisition

Data analysis

Instrument shutdown

Instrument Startup

Startup function:

Warms up the lasers

Initializes the pumps

Primes the instrument fluidics

Flushes out the shutdown solution

Ensures that:

All fluidic lines are clean

The fluidic lines and the syringe pump are filled with fresh focusing fluid

The lasers are warmed to operating temperature

Startup



Takes 3.5 minutes

Uses ~25 ml of focusing fluid

Instrument Startup

From power off

1. Turn on autosampler remove cleaning plate and close tray door
2. Turn on the cytometer and wait until status light remain a solid, bright blue
3. Open the Attune NxT Software
4. Open the instrument tab and click startup icon

From Sleep mode

After shutdown completes

1. Close/re-open the software
2. Login
3. Open the instrument tab
4. Open the instrument tab and click startup icon

Once Startup is complete, the status lights are solid green and the status bar displays the *Ready* icon



Instrument Start-up

Performance test

Create experiment

Experiment settings optimization

Compensation

Data acquisition

Data analysis

Instrument shutdown

Performance Test

- Allows you to monitor performance of the instrument
- Critical to ensure accuracy and sensitivity of instrument
- Provides information about the lasers and detection channels
- Run immediately after launching the software application and running **Startup**

There are 2 parts to Instrument Performance Tracking:

1. Baseline Calculation (BL)
2. Daily Performance test (PT)

Performance Test 'when to run'

Daily performance Test

Run daily - everyday samples are run/recorded

Baseline Calculation (advanced user/admin/sysadmin only)

Performed at time of installation by Field Service Engineer (FSE)

After any major service (FSE)

Every time the bead lot changes (User)

When recommended by FSE or FAS

Attune® Performance Tracking Beads

- A mixture of beads of four fluorescence emission intensities in equal concentration


- Blank
- Dim
- Medium
- Bright

- 3mL vial

PN: 4449754

Lot # 2029773xx

Lot # batch #



3 drop of beads per 2 ml of Focusing Fluid or PBS

Notes: Data for a new lot of beads can be downloaded from the Performance Tracking beads webpage on the ThermoFisher website


Documents

Product literature

- [Attune \(not NxT\) Performance Tracking Bead Lot File--Lot 1759476](#)
- [Attune \(not NxT\) Performance Tracking Bead Lot File--Lot 756080](#)
- [Attune NxT Performance Tracking Bead Lot File Installer--Lot 1759476](#)
- 📄 [Brochures & Specifications: Product Information Sheet: Attune™ performance tracking beads](#)
- 📄 [Detecting Human Circulating Endothelial Cells Using the Attune® Acoustic Focusing Cytometer](#)

Performance tracking beads

Bead doublets and multiples



Size:

small

large

large

large

Dye level:


none

low

med

high

Performance tracking beads



Size:

small

large

large

large

Dye level:

none

low

med

high

Baseline Calculation and Performance Test

- **Baseline** uses Performance Tracking beads specifications (MESF values) to define **initial** status of the Attune® NxT™ Cytometer
 - PMT voltages are adjusted to place the brightest bead at target MFI values; voltage value for each channel is reported
 - The robust % coefficient of variation (Robust %CV) of the brightest bead is recorded
 - Relative quantum efficiency (Qr) of each detector is determined
 - Relative Background level (Background) of each detector is determined
 - Linear regression (Linearity) is calculated and recorded
 - Area scaling factor (ASF) is calculated and reported for every laser and automatically updated in “Advanced Settings”
 - Laser delay setting is automatically calculated
- Performance test uses Performance Tracking beads to monitor changes over time
 - Same process/measurements as Baseline - reports all
 - The change of the PMT voltage (**Δ PMT**) from baseline is also reported

Baseline Report

✓ Baseline test successful

Baseline 756080D - 7/24/2014 < 7/24/2014 12:03:42 PM

Channel	PMTV	Target MFI	MFI	Robust %CV	Qr	Background	Linearity	ASF	Laser Delay	Result
FSC	575	300000	302004	1.28 %	0.000	0	0.000	1.02	1100	✓
SSC	358	300000	306486	3.03 %	0.000	0	0.000	1.02	1100	✓
BL1	381	300000	300836	1.30 %	0.060	101	1.000	1.02	1100	✓
BL2	361	300000	304638	1.57 %	0.058	135	0.965	1.02	1100	✓
BL3	413	300000	301436	1.84 %	0.051	37	1.000	1.02	1100	✓
RL1	367	300000	315794	3.99 %	0.064	40	0.998	0.97	1557	✓
RL2	378	300000	309933	3.72 %	0.013	176	1.000	0.97	1557	✓
RL3	407	300000	310010	3.59 %	0.079	77	0.997	0.97	1557	✓
VL1	297	300000	291418	0.97 %	0.014	949	1.000	0.81	698	✓
VL2	385	300000	304298	1.11 %	0.021	224	0.998	0.81	698	✓
VL3	375	300000	303931	1.32 %	0.023	98	0.996	0.81	698	✓
VL4	433	300000	312414	2.21 %	0.006	235	0.984	0.81	698	✓
YL1	401	300000	301724	1.90 %	0.110	40	0.999	0.71	239	✓
YL2	390	300000	308583	1.71 %	0.071	38	0.973	0.71	239	✓
YL3	430	300000	300002	2.12 %	0.030	100	0.999	0.71	239	✓
YL4	501	300000	302131	3.08 %	0.004	320	1.000	0.71	239	✓

Pass ✓

Fail ⚠



Performance Test Report

Baseline		756080D - 7/9/2014								7/23/2014 3:27:05 PM	
Channel	PMTV	Delta PMTV	Target MFI	MFI	Robust %CV	Qr	Background	Linearity	ASF	Laser Delay	Result
FSC	568	-9	300000	299756	2.06 %	0.000	0	0.000	1.02	1100	✓
SSC	350	-10	300000	288385	3.89 %	0.000	0	0.000	1.02	1100	✓
BL1	379	-1	300000	300519	1.24 %	0.060	115	1.000	1.02	1100	✓
BL2	359	1	300000	306803	1.34 %	0.057	118	0.965	1.02	1100	✓
BL3	410	-3	300000	297859	1.99 %	0.052	45	1.000	1.02	1100	✓
RL1	366	-2	300000	319156	3.83 %	0.070	42	0.998	0.96	1563	✓
RL2	374	-6	300000	292375	3.76 %	0.012	159	1.000	0.96	1563	✓
RL3	407	-3	300000	303972	3.82 %	0.059	72	0.997	0.96	1563	✓
VL1	301	4	300000	297676	1.58 %	0.008	579	1.000	0.82	694	✓
VL2	383	-3	300000	297463	1.15 %	0.020	282	0.998	0.82	694	✓
VL3	374	-6	300000	300917	1.40 %	0.023	93	0.995	0.82	694	✓
VL4	429	-8	300000	289260	2.17 %	0.005	221	0.984	0.82	694	✓
YL1	400	-3	300000	292582	1.44 %	0.092	34	0.999	0.68	229	✓
YL2	390	-3	300000	304246	1.31 %	0.067	36	0.973	0.68	229	✓
YL3	430	-3	300000	294263	2.16 %	0.028	94	0.999	0.68	229	✓
YL4	500	-3	300000	288991	3.16 %	0.005	309	1.000	0.68	229	✓


Pass
















Fail



Performance Test Territory

 Performance test completed with errors

Baseline: 2029773 - 1/22/2020 2/25/2020 8:36:42 AM

Channel	PMTV	Delta PMTV	Target MFI	MFI	Robust %CV	Qr	Background	Linearity	ASF	Laser Delay	Result
FSC	0	0	300,000	0	0.00 %	0.000	0	0.000	1.20	1100	
SSC	339	-3	300,000	309,302	4.17 %	0.000	0	0.000	1.20	1100	
BL1	444	-1	300,000	305,044	1.77 %	0.053	132	1.000	1.20	1100	
BL2	408	0	300,000	303,814	1.68 %	0.055	193	1.000	1.20	1100	
BL3	439	-3	300,000	301,628	2.37 %	0.057	36	1.000	1.20	1100	
VL1	319	12	300,000	300,410	0.91 %	0.034	1724	1.000	0.83	732	
VL2	348	13	300,000	305,712	1.06 %	0.030	448	0.997	0.83	732	
VL3	400	15	300,000	305,947	1.14 %	0.039	67	1.000	0.83	732	
VL4	486	16	300,000	299,237	2.22 %	0.006	201	0.990	0.83	732	
YL1	423	-1	300,000	300,913	1.73 %	0.129	104	0.999	1.00	368	
YL2	366	0	300,000	304,590	1.94 %	0.106	42	1.000	1.00	368	
YL3	421	-1	300,000	298,878	2.71 %	0.030	171	1.000	1.00	368	
YL4	510	-1	300,000	300,466	3.17 %	0.006	240	1.000	1.00	368	

Pass



Fail





Pass - All statistics and calculations for the channel meet the criteria set by the Baseline calculation.



Fail - One or more of the statistics or calculations for the channel deviate significantly from the target set by the Baseline calculation.


e.g. Δ PMT exceeds 100 mV

%rCV – detector specific but ranges from 3-5%

Levey-Jennings Report


- Tracks “Performance Test” data
- Provides a visual indication of the cytometer performance over time.
- Monitor shifts and trends in cytometer performance and

Select: statistic



Select laser to view


Scatter Channels:



View/Print individual reports

Performance History Report

- Pass/fail status of all PT runs
- Filter date range
- Quick visualization of 'health' of instrument
- Click on result to open associated PT test
- Max 180 days for current baseline
- Not available if no baseline exists



Instrument Start-up

Performance test

Create experiment


Experiment settings optimization

Compensation

Data acquisition

Data analysis

Instrument shutdown



Create a new experiment by selecting New Experiment

New Experiment

Experiment type: Experiment

Tube Experiment(2)

Tube

Plate

Use workspace:

Load Default workspace

Use instrument settings:

Load Default instrument settings

Create 1 group(s) for this experiment

Create 1 tube samples for each group

Notes:

OK Cancel

Select Experiment type and click OK

Instrument Start-up

Performance test

Create experiment

Experiment settings optimization


Compensation

Data acquisition


Data analysis

Instrument shutdown

Experiment settings optimization





Instrument
settings panel



1. Parameters

- Expand parameters section
- Optional: Select customized configuration
- Deselect detectors/channels not required
- Select A – H – W as needed
- Select **label** from drop down menu
- Enter target name
- After information has been entered, collapse the section



Create Workspace

The screenshot displays the Attune™ NxT Software v2.5 interface. The main window is titled "Preview - All Events" and shows a scatter plot of SSC-A (10³) versus FSC-A (10³). The plot is labeled "All Events - Sample". The left panel shows the "Collection Panel" with a circular progress indicator at 0%, "0 events", and "Sample remaining". Below this are buttons for "Stop", "Start Up", and "Clear". The "Run Protocol" section includes "Apply to experiment", "Set as default", and "Load..." buttons. The "Flow Options" section shows "Acquisition Vol" set to 50 µL (100 µL Total Draw Volume) and a flow rate of 200 µL/min. The "Stop Options" section shows a threshold of 10,000 events on "All Events". The right panel shows the "Instrument Settings" section with a table of parameters.

Instrument Settings

Enabled	Target	Label	A	H
<input checked="" type="checkbox"/>	FSC		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	SSC		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	BL1		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	BL2		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	BL3		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	RL1		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	RL2		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	RL3		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	VL1		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	VL2		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	VL3		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	VL4		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	YL1		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	YL2		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	YL3		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	YL4		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>



The bottom status bar shows "User Time: 00:10", "Workspace: ExperimentT1", "Plate Time:", "Instrument Settings: Experiment", and a zoom level of 100%.

2. Adjust PMT Voltages

- Expand the **Voltage** section
- Adjust FSC & SSC voltages to position cell population on the scatter plot
- Adjust Fluorescence Channels voltages based upon the positive control
- If you do not have a positive control, set the voltage for the negative population so that the events are displayed between 100-1000 relative fluorescent units


Note:

- It is recommended to perform a voltage walk and calculate the staining index for each fluorophore
- If compensation controls are recorded, all the fluorescence channels voltages are disabled (i.e. grayed out)




Match the FSC/SSC voltage settings with the plots

LG1




FSC 400
SSC 400

PT beads




FSC 450
SSC 350

PT beads



FSC 350
SSC 350

PT beads





FSC 500
SSC 325

3. Adjust Threshold (optional)

- Expand the **Threshold** section
- Way to get rid of unwanted events (i.e. noise) before sample has been recorded
- Default setting: OR FSC 25 x 1000.
- Can be set on a single or up to 4 scatter and/or fluorescence channels
- Data not meeting threshold criteria is permanently lost

- Exclude coincident events option:
 - checkbox
 - # of aborted events in FCS file header


At least 1 threshold must be set . If all are set to OFF, no data is displayed on the workspace




3. Adjust Threshold (optional)

A threshold: is a way to get rid of unwanted events (e.g. debris) before a sample is recorded


**Above threshold:
an “event”**



**Below threshold:
not an “event”**



Examples of Threshold Adjustment



Stop Option – 10,000 events on All Events

Advanced instrument settings

***(Administrator and Advanced User
permission required)***

- Width threshold setting
- Area scaling factor (ASF)
- Front + Rear window extensions

The screenshot shows the 'Advanced Settings' window. At the top, there is a checkbox labeled 'Automatically update Experiment level IS' which is checked. Below it are two buttons: 'Set as Experiment Level IS' and 'Set as Default'. A dropdown menu is open, showing a list of settings: 'Parameters', 'Voltage', 'Threshold', and 'Advanced settings' (which is highlighted). Below the dropdown, the 'Width threshold setting (x1000)' is shown with a slider and a text box containing '1.0'. The 'Area scaling factor' section has four rows: 'Blue' with a value of 1.07, 'Red' with 1.08, 'Violet' with 1.14, and 'Yellow' with 0.99. The 'Window extension setting' section has a row for 'Front + Rear' with a value of 0.


Color	Area scaling factor
Blue	1.07
Red	1.08
Violet	1.14
Yellow	0.99

Setting	Value
Width threshold setting (x1000)	1.0
Window extension setting - Front + Rear	0

More information can be found in the
SW User Guide Rev C (loaded on the desktop)

Area Scaling Factors

- The Area Scaling Factor (ASF) is a correction parameter that sets the area and height measurements at parity:



- Area Scaling Factors are Calculated and Applied from the results of the performance test (V2.4 – not V2.2)

Channel	PMTV	Target MFI	MFI	Robust %CV	Qr	Background	Linearity	ASF	Laser Delay	Result
FSC	555	300000	302000	1.80 %	0.000	0	0.000	1.01	1100	✓
SSC	318	300000	304325	3.44 %	0.000	0	0.000	1.01	1100	✓
BL1	403	300000	301618	1.34 %	0.051	86	1.000	1.01	1100	✓
BL2	348	300000	311511	1.34 %	0.068	138	0.966	1.01	1100	✓
BL3	385	300000	305088	1.72 %	0.047	18	1.000	1.01	1100	✓
RL1	480	300000	304847	7.33 %	0.016	59	0.994	0.78	1512	✓
RL2	493	300000	307636	7.70 %	0.002	448	0.999	0.78	1512	✓
RL3	510	300000	307720	7.87 %	0.007	217	1.000	0.78	1512	✓
V/L 1	306	300000	304797	1.07 %	0.015	2479	1.000	1.01	676	✓

Area scaling factor

Blue 1.2


Red 1.1

Violet 1.2

Yellow 1.0

Area Scaling Factors - Why would I change them?

ASF changes as a function of particle size




Adjusting ASF

If Area (MFI) is greater than Height (MFI); enter a value <1









If Area (MFI) is less than Height (MFI); enter a value >1

*Note: Adjustment must be made **prior** to recording data*








Examples of Adjusting Area Scaling – 1st Method

Compare FSC-H vs FSC-A in histograms



Examples of Adjusting Area Scaling – 1st Method

Compare FSC-H vs FSC-A a bivariate plot



1. Draw a diagonal line
2. Cells should fall as close to the line as possible
3. But towards the right of the line
4. AS 1.0 and 1.1 are optimal for this sample
5. Cell lines usually need to be adjusted

Instrument Start-up

Performance test

Create experiment

Experiment settings optimization

Compensation


Data acquisition

Data analysis

Instrument shutdown

Why do we need to compensate?

- Every fluorescent molecule emits light with a particular spectrum unique to that molecule
- These emission spectra overlap and in some cases is very significant
- *Compensation* is the process by which we correct for "spillover"




<https://www.thermofisher.com/order/spectra-viewer>


How does fluorescence 'spillover' look in data plots?

Single stained sample: **PE-Cy5.5**

Uncompensated




Compensated



Compensation is the process of correcting the **spillover** from our signal (e.g. PE-Cy5.5) into YL4 and each secondary channel into which it is detected/measured.

How? We match median fluorescence of positive and negative populations

Compensation Confirmation



Parameters: YL4-H vs YL1-H

Gate: CELLS

Experiment: 14H13_INKT

Group: INKT MATRIX


Sample: NOT-IN_STATS

Name	Count	%Gated	%Total	Y Median
■ All Events	1106537	100.00	100.00	413.88
■ -/APC-CD8+(1)	1	0.00	0.00	4104.51
■ +/-APC-CD8+(1)	2	0.00	0.00	3967.83
■ -/APC-CD8-(1)	515380	59.60	46.58	419.63
■ +/-APC-CD8-(1)	349313	40.40	31.57	410.26

- Check Median of fluorescence is equivalent in non-targeted channel
- Check for all dual parameters plot combination

Compensation and Spreading Error

PE into PerCPCy5.5



Basic Rules of Compensation

1. Unstained cells

- Only required if sample does not have an unstained population
- Background fluorescence should be the same for the positive and negative populations

2. Single color controls for each fluorochrome

- Compensation color must be matched to your experimental color (FITC cannot substitute for GFP)
- The actual tandem dye being used in the sample staining must be used in the single-color control. Lot numbers must match.

3. Controls need to be at least as bright as the brightest positive sample


- Titrate antibody on compensation beads. Beads are efficient at antibody binding and maybe too bright and off scale compared to cell staining
- Log separation between the negative and positive peaks
- Brightest in the primary detector (its own channel) or spilling over into other channels less than 100%

4. Collect enough events.

Compensation Set up


In *Compensation* tab, click on **Compensation Setup**. Or

Double-click on Compensation in Experiment Explorer



Compensation setup options:

- Source: Tube
- Parameter: Area or Height
- Autofluorescence:
 - Negative gate,
 - Unstained control
 - None
- Fluorescent channels



Compensation: Background Fluorescence Modes


Auto Fluorescence Correction Choices:

- Negative Gate
- Unstained Control
- None

Background Mode	When to Use?
Negative Gate	With different controls such as cells and beads; or using different cell populations (lymphs and monos).
Unstained Control	When all controls are of the same type (beads, all lymphs)
None	Rarely used but in cases where background is negligible or cannot be ascertained.

Spillover Matrix

- At the end of Auto-compensation:
 - Spillover Matrix is automatically calculated
 - Compensation is applied to all samples



Spillover						
	BL1-H	BL3-H	RL1-H	RL3-H	YL1-H	YL4-H
BL1-H	100.00	0.31	0.12	0.14	0.04	0.04
BL3-H	1.76	100.00	1.65	0.67	81.71	14.70
RL1-H	0.10	0.29	100.00	22.10	0.00	1.08
RL3-H	0.17	0.00	0.39	100.00	0.01	3.64
YL1-H	1.29	3.68	0.09	0.02	100.00	0.50
YL4-H	3.24	0.22	0.11	34.20	2.44	100.00

Spillover Matrix

The spillover matrix will assign a numeric value of percent spillover of a fluorophore into other detectors once all controls have been recorded. To read the matrix, consider the column on the left the fluorophore (ie FITC) and the top row the detector.

If FITC is the fluorophore assigned to BL-1 detector, then the table below shows that FITC spills into the BL-1 detector 100%, the BL-2 detector 38.10% and the BL-3 detector 5.86%


Spillover	BL1-A	BL2-A	BL3-A	RL1-A
BL1-A				
BL2-A				
BL3-A				
RL1-A				

	BL1-A	BL2-A	BL3-A	RL1-A
BL1-A	100.00	38.10	5.86	0.01
BL2-A	0.51	100.00	24.18	0.00
BL3-A	0.22	25.40	100.00	13.92
RL1-A	0.06	0.04	0.01	100.00

On Plot Compensation Adjustment Tools

NOTE: Be very careful when using manual compensation as it may poorly effect your results.



Tool allows manual adjustment of a selected plot



% Spillover

Quadrant statistics

Dragging populations



Instrument Start-up

Performance test

Create experiment

Experiment settings optimization

Compensation

Data acquisition


Data analysis

Instrument shutdown

Immunophenotyping Example: Three Color Experiment

Sample: Leukocytes from Human Blood




Measure: % T-lymphocytes CD4+
% T-lymphocytes CD8+






	Antibody	Fluorescent Probe
1	Anti-CD3	Alexa Fluor® 488
2	Anti-CD4	R-PE
3	Anti-CD8	R-PE Alexa Fluor® 700 dye tandem


→ T-lymphocytes specific

Three Color Experiment – Data presentation




Three Color Experiment – Data Representation





Instrument Shutdown




Shutdown Function

- Sanitizes the instrument
- Cleans and rinses the fluid lines
- Refills fluid lines with shutdown solution
- Requires 10% Bleach freshly prepared

Ensures that:

- Fluid lines refilled with a solution that prevents crystal formation and bubbles.


Instrument Shutdown



- In *Instrument* tab, select **Shutdown** function

- Select:

Quick	15 cycles / 30 minutes
Standard	25 cycles / 45 minutes
Thorough	35 cycles / 60 minutes



- Follow the on screen instructions
- Once *Shutdown* starts, *log out of the software if necessary*
- Use 10% bleach - freshly prepared
- At the end of the Shutdown script:
 - the Cytometer will be in 'dream state'
 - the Autosampler will be in standby.
 - The fluid in the tube on the SIP will be shut down solution, not bleach solution

Attune® NxT Cytometer System Maintenance


Instrument cleaning guide

Between samples	<ul style="list-style-type: none">• Rinse – automatically initiated when SIP is lowered (for tubes), or set in <i>run protocol</i> for plates• Sanitize SIP between sticky samples or cell counts
Between users / experiments Use: 1) if there is ≥ 30 min between users. 2) If there is < 30 min between users.	<p>1) Unclog then Quick Deep Clean - 30 minute cleaning routine (click on the arrow below the Deep Clean icon to select <u>Quick</u>)</p> <p>or</p> <p>2) Unclog then 2X Sanitize SIP / Sanitize Autosampler SIP (plate experiments) –</p> <ul style="list-style-type: none">1st time with 3 mL 10% Bleach2nd time with 3 ml Wash or De-bubble solutions
End of day (3 steps)	<ul style="list-style-type: none">• <i>Unclog</i>• <i>**SIP Sanitize with 1:3 dilution of Attune Cleaning solution</i>• <i>Thorough Shutdown</i> (click on the arrow below the Shutdown icon to select <i>Thorough</i>)

Note: Always wipe the outside of the SIP after doing a SIP Sanitize



Daily - Visual inspection

- Fluidics compartment: make sure there are no fluids or salt residues on the floor of the compartment, around the connectors, or tube junctions
- Check the fluids level. Fill/empty as needed:
 - Focusing fluid
 - Wash solution
 - Shutdown solution
 - Waste
- Visually inspect the SIP




Daily - Visual inspection

- **Syringe compartment** – confirm there is no fluid or salt residue on the floor of the compartment
- **Syringe** – finger tighten the syringe; change if there is a leak or salt residue builds-up
- **Focusing fluid filters** – located behind the wash and shutdown fluidics bottles. Change if there are any signs of debris/dirt or if the sample pump stays on too long



Filling (or emptying) Fluid Tanks



1. Remove the sensor cable from the instrument
2. Press the metal release buttons to free the tubing
3. Fill or empty as needed with RT solutions
Large tanks – 1.8 L
Small tanks – 175 mL
4. Return to cytometer and reconnect the fluid line then the sensor line

	Fluid line	Sensor cable
Disconnect:	2	1
Reconnect:	1	2

IMPORTANT

1) Connecting the sensor cable while leaving the fluid line disconnected may result in increased back pressure and introduction of air into the system.



2) The Attune® NxT Acoustic Focusing Cytometer must be idle before refilling the fluidics containers.

3) Do not pull on the lines.




Potential problem – Contamination

- Check tanks for cloudiness or debris in the solutions or brown marks on the sensor
- Fill the emptied waste container with a **full strength bleach up to the bleach fill mark (bottom line) on the bottle



Fluidics Replacement Parts



- Replacement part's:
 - 1.9 L waste tank # 100022156
 - 1.9L focusing fluid solution tank # 100022155
 - 175 ml wash solution tank # 100022151
 - 175 ml shutdown solution tank # 100022154
 - AAS focusing fluid solution tank # 4477847
 - AAS waste tank # 4477850

System Decontamination



- Function which facilitates the decontamination of the Attune® NxT™ Acoustic Focusing Cytometer and the Attune® Autosampler fluidics.
- Sanitizes the system and fluidics bottles with Bleach and Wash solutions.
- Mostly automatic - 60 minutes, 3 phase operation with on screen instructions
- **REPLACE FOCUSING FLUID FILTERS**


Note: this function is only available to advanced users and administrators

When?

- As a **quarterly** maintenance routine to prevent and reduce microbial growth within the instrument and fluidic bottles
- If the system is likely to be idle for more than two weeks (run it in place of the Shutdown function)
- If the instrument has been idle for more than two months
- Anytime contamination in the fluid lines is suspected – i.e. event rate is too high
- Prior to any service work or shipment for service



NEW: Improved instructions for System Decontamination

- Decontaminate System is divided in to three phases: **Sanitize, Rinse and Refill**



This screen is displayed when bottles are disconnected

Phase 1



Decontaminate System

Phase 2

Decontamination - Phase 2 of 3

Sanitize	Rinse	Refill
----------	-------	--------

Rinse out all fluid bottles with deionized water
Fill Attune™ and autosampler focusing fluid bottles with 500mLs of deionized water
Fill shutdown bottle with 100mLs of deionized water
Fill wash bottle with 100mLs of deionized water

Previous Next Cancel

Decontamination - Phase 2 of 3

Sanitize	Rinse	Refill
----------	-------	--------

Ensure all fluid lines and cables are reconnected
Load a clean, empty tube and raise tube lifter

Previous Next Cancel

Phase 3

Decontamination - Phase 3 of 3

Sanitize	Rinse	Refill
----------	-------	--------

Lower tube lifter
Remove plate from autosampler

Previous Next Cancel

Decontamination - Phase 3 of 3

Sanitize	Rinse	Refill
----------	-------	--------

Rinse out all fluid bottles with deionized water
Replace all fluids in all bottles with appropriate solutions

Previous Next Cancel

Decontamination - Phase 3 of 3

Sanitize	Rinse	Refill
----------	-------	--------

Reconnect all fluid lines and bottle cables
* Replace both filters (located behind the fluid bottles)

Previous Next Cancel




After System Decontamination: Prime Fluidics

After System Decontamination and replacement of focusing fluid filters, the Attune NxT must be primed to remove air from the new filters

To Prime the Fluidics System:

1. Run 3 Startup cycles
2. Run 2 Debubble cycles using Debubble solution or wash solution
3. Run 1 Rinse cycle
4. Run Performance Test. If PT fails, run two debubble cycles and repeat Performance Test.

Your system is now ready to use



Maintenance – Focusing Fluid Filters







Focusing fluid filter
PN 100022587



- Two filters located behind the Wash and Shutdown bottles
- The filters may grow some contamination over time. If discoloration is evident, replace the filters.
- Replacing focusing fluid filters **quarterly** reduces the risk of any potential contamination in the lines. **ALWAYS REPLACE AFTER SYSTEM DECONTAMINATION.**

How to change Focusing Fluid Filters?

- Remove shutdown and wash bottles from bottle bay
- Unscrew the top luer fitting
- Unscrew bottom luer fitting and remove filter
- Put fittings on new replacement filter and re-attach to unit
- Ensure arrow on filter points in the direction of fluid flow (down)
- **Prime the filters by running 3X Startups and 2X Debubble**



Maintenance – Sample Syringe



- Potential problems:
 - Check for leaks
 - Erratic or no fluid draws up from SIP
 - Erratic or no fluid draws up from fluidics tanks
- Part numbers:
 - 1 ml syringe # 100022591
- Replacing Frequency: as needed, at least bi-annually

How to change the sample syringe?

- Run Shutdown
- Open the Syringe Pump door located on the left side of the cytometer
- Loosen the knurled thumbscrew below ball end of syringe
- Unscrew top portion syringe from valve head
- Remove syringe from ball end, pull out and replace with new syringe
- Tighten the syringe 1/4 turn past initial contact of the Teflon insert into the valve for a liquid seal.
- Properly seat the ball end of the syringe. Tighten the knurled thumbscrew below the ball end.

Notes: No tools should be used to tighten the syringe to the valve



What does the Attune NxT Database utility do?

- 1. Backup User Data:** backups the entire database, files plus database data to a folder a single time (duplicates everything upon initiation)
- 2. Restore User Data:** restores the contents of the folder from a backup so that the backup becomes the current version of the database (eg, the “live” database)
- 3. Automated backup:** schedules automatic backup of user data and the database
- 4. Re-Install Database:** used to reset the database to the **new, no data added state**. This should **rarely need to be used**
 - When would this be suitable? If recommended by FSE
 - **Please reference detailed slide deck for more information**

Administrator

Set up user accounts with operator privileges
User, Advanced User, Administrator

Back up experiments to secondary storage

Virus protection – scan thumb drives before connecting to Attune computer

Defragment the hard drive *weekly*

Network connection (optional)

Maintaining computer efficiency & data quality

- **Operator/user**

De-select those parameters not needed - minimize file size

***Do not clutter the Experiment Explorer ***

Close experiments not currently active (collapse all)

Export then delete experiments from the browser

Export experiments - .atx or .apx

FCS files – 3.0 or 3.1 format

If you need to scroll through the list of experiments in the experiment explorer – it is time to backup and delete

Virus protection – scan thumb drives before connecting to Attune computer

Uninterrupted Power Supplies

We recommend the use of a 1.5-kVA uninterruptible power supply (UPS), especially in areas prone to Power failure.

Anti-Virus software –

- Disable or deactivate antivirus software and antispyware during use of the Attune® NxT Acoustic Focusing Cytometer.
- Antivirus and antispyware monitoring can interfere with data collection, resulting in data loss

Maintenance Summary

Procedure		Frequency
Startup and Shutdown		Daily (2-3 x per week)
Visual inspection of	sample injection port (SIP), fluidics tanks and connections syringe pumps	Daily
Fluidics maintenance – cleaning routines on Instrument tab <i>see the posted daily cleaning guidelines</i>	Rinse	Daily – between samples
	Sanitize SIP	Daily – between samples, experiments, or users
	Deep Clean	Daily – between experiments
	Debubble	Daily – as needed
	Unclog	Daily – as needed
Sanitize SIP with Attune Cleaning solution	Usage - heavy: ≥6 hr per day, ≥6 plates per day, multi-lab shared system	Daily before shutdown
	Usage - light to moderate <6 hr per day	Once a week before shutdown
Fluid bottles	Empty, clean, rinse, refill	Monthly
System Decontamination		Every 3-6 months
Replace focusing fluid filters		Every 3-6 months
Syringe replacement		6 months or as needed
Optical filter and mirror inspection		Monthly or less
Secondary backup of experiments	Export experiments (.atx or .apx) or FCS files	Daily, weekly
Attune Database	Schedule auto backup	Daily, weekly or monthly

Service Preventive Maintenance (PM)

- Replacement of:
 - Focusing Fluid Filters
 - 1 mL Syringe
 - Sample probe (SIP)
- Cleaning of:
 - Interior and Exterior of unit
 - Optical filters
- Check for leaks
- Check fittings and valves
- Run a system diagnostics test
- Check pinhole and laser alignment
- Computer Maintenance

Ensures maximal performance



√ Included ○ Option	AB® Complete	AB® Assurance	AB® Maintenance
Planned maintenance	√	√*	√
On-site service--Labor	√	√	
On-site service--Parts	√	√	
On-site service--Travel	√	√	
Remote instrument monitoring diagnostics	√	√	
Telephone Support (within 3 hours)	√	√	√
Application technical support	√	√	
On-site application consulting	√		
Qualification service	√	○	○
Computer System Validation	○	○	○
On-site response time	Guaranteed next business day	Typical 2 business days	Guaranteed 3 business days**

A service plan from ThermoFisher Scientific can help you:

- Maximize productivity
- Optimize your laboratory's efficiency
- Lower the cost of ownership
- Obtain unmatched availability of critical laboratory systems
- Increase quality
- Lower costs by minimizing lost data, samples, or reagents

* Available with 1 or 2 pm/year

** After purchase order has been received by ThermoFisher Scientific

Changing the location of the Attune® NxT?

Bench space needed: W x H x D

Width: 50 inches

- 22.9 in. for Attune® NxT

- 3.75 in. for access to side syringe compartment

- 22.5 in. for Computer system

Height: 29 inches to allow the hinged lid to fully open

Depth: 23.1 inches

- 17.1 in. for Attune® NxT cytometer

- 2.5 in. for adequate ventilation behind the instrument

- 4 in. for the fluidics bottles in front of the unit

Logins - all case sensitive

Instrument is powered on in a set order:

- 1) NxT Autosampler
- 2) NxT Cytometer
- 3) Open Attune NxT software

Instrument:

User name: INSTR-ADMIN

Password: INSTR-ADMIN

Attune NxT Software:

User name: admin

Password:

Individual users (account can be setup by someone with administrator privileges):

User name: avoid ALL CAPS

Password: same as user name. Temporary until first login

Thank you!

Please provide feedback if/when surveyed

©2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

*The following are reference slides for Attune NxT software

- Hardware User Guide



https://assets.thermofisher.com/TFS-Assets/LSG/manuals/100024235_AttuneNxT_HW_UG.pdf

- Software User Guide



https://assets.thermofisher.com/TFS-Assets/LSG/manuals/100024236_Attune_NxT_SW_UG.pdf

- Maintenance and Troubleshooting Guide



<https://assets.thermofisher.com/TFS-Assets/BID/Product-Guides/attune-nxt-maintenance-troubleshooting-guide-quick-reference.pdf>




Panels



Ribbons and Tab (#1)




Home



View




Instrument




Autosampler

More information in the SW user guide


Ribbons and Tab (#2)




Workspace



Compensation





Overlay builder



Statistics

More information in the SW user guide

Setting User Preferences



Colors and Themes

Gate and Sample Overlays

Group Color

Experiment Color

Export Options

Plot Export Options

Default Plot Image Format

Export Size (scaled to width)

Statistics Options

Header

Statistics

Concentration Units

Stats Box Style

Decimal Settings

Fonts and Styles

Plot Title and Axis Label Style

Text Box Style

Gate Options

Gate Style

Gate Label Style

Gate Label Naming Options

Quick Select Gate Names

Workspace

plots

gate

scales - linear, log and hyperlog

statistics box and results table

Customize

Heat Map View and Setup

Filter Configuration Manager

Workspace - Plots







Histogram – single parameter plot showing number and distribution of events

Dot plot – two parameters plot where each axis represents the signal intensity of one parameter

Density plots – two parameters plot where colors represent the density of a population of events with the same intensity

Density precedence – a combination of Dot and Density display. A gradient is used to indicate the number of events within each of the plot bins and color is used to display the parent gate of events present.



Use the *Customize* Panel to change plot titles, axis labels, axis scaling, plot type etc.

Note: *Defining the workspace and gating strategy prior to data acquisition is highly recommended. Final gate adjustment should be made after a file is recorded.*

Workspace – Gating Tools




- Regions and gates are commonly used in data analysis to identify population subsets



- Gate** is a shape or object that is drawn around a population of interest on one or two parameter plots
- Region** is defined when gates are used to isolate a specific group of cytometric events from a large set of data
- Gates are displayed in a hierarchy or family tree
- Region may be exported as an fcs file by right-clicking on the gate and choose Export to fcs file

Daughter Plots


- Displays data only from the region selected



- To create
 - Right click on the gate and select **Daughter plot** and the type of plot
 - OR
 - Right click on the plot and select **Set population**
 - OR
 - Left click on the Hyperlink in the plot title and select gate name from the drop down list
 - Note: do not use quadrant gate as a parent gate

Workspace - Derived gates

- Gates can be customized by using Boolean logic (OR, AND, NOT, XOR) to link multiple gates together



The screenshot shows the 'Gate Logic' dialog box. It has a title bar with a close button. Inside, there are four fields: 'Gate Name:' with the text 'CD8- Lymphocytes'; 'Gate Color:' with a black color swatch; 'Gate Definiton:' with a dropdown menu showing 'Lymphocytes'; and another dropdown menu showing 'AND NOT' followed by a text field containing 'CD8+'. At the bottom left is a button with three dots, and at the bottom right are 'OK' and 'Cancel' buttons.

AND gates = all events that are shared

OR gates = all events found within 2 or more individual gates

NOT gates = all events found outside the gate

XOR gates = unique events found within an individual gates


Note: When naming a derived gate with two words, use parentheses to enclose both words.

Derived gates can only be created using regions.

Impact the way data is visualized

- Linear – data spread over a single order of magnitude
 - FSC
 - SSC
 - DNA content
- Log – data spread over a wide range (>1 order of magnitude)
 - Fluorescent channels
 - FSC +/- SSC – bacteria, small particles and blood samples


- **Log scale** – cannot correctly represent values for cells whose fluorescence values fall at/or below zero - compensated data piles up on the axes



Negative fluorescence results from

- Background subtraction
- Compensation

- HyperLog™ – similar to LinLog scale used in classic Attune®
 - Logarithmic scale at the high end
 - Transitions to linear scale in the region around zero




- Use to determine if compensation is correct
 - Correct – double neg population will distribute symetrically around autofluorescence level
 - Overcompensation – double negative population center below autofluorescence value
 - Undercompensation – double negative population centers the distribution above autofluorescence value


Modern Flow Cytometry: A Practical Approach Clin Lab Med. 2007 September ; 27(3): 453–v.

Visual

Workspace plots



Heat Map



Statistics

Statistics box (on plot)

Name	Count	%Total	%Gated	Concentration
▢ All Events	25,218	100.000	100.000	2,802.000
▢ Cells	20,002	79.316	79.316	2,222.444
▢ Live	18,263	72.420	91.306	2,029.222
▢ RBC	12,163	48.231	66.599	1,351.444
▢ live wbc	5,435	21.552	29.760	603.889
▢ CD45R	1,018	4.037	18.730	113.111
▢ CD11b	3,665	14.533	67.433	407.222


Results table

Filename	Count	Concentration	%Total	%Gated
test000.fcs	22801	1341.24	100.00	100.00
test000.fcs	22801	1341.24	100.00	100.00
test000.fcs	22801	1341.24	100.00	100.00
test000.fcs	22801	1341.24	100.00	100.00
test000.fcs	22801	1341.24	100.00	100.00
test000.fcs	22801	1341.24	100.00	100.00
Gate: CD11b				
cont1.fcs	3665	407.22	14.53	67.43
test1.fcs	4147	243.94	17.94	72.91
test000.fcs	3482	204.82	15.27	73.27
Gate: CD45R				
cont1.fcs	1018	113.11	4.04	18.73
test1.fcs	395	23.24	1.71	6.94
test000.fcs	306	18.00	1.34	6.44
Gate: Cells				
cont1.fcs	20002	2222.44	79.32	79.32
cont1.fcs	20002	2222.44	79.32	79.32
test1.fcs	20002	1176.59	86.55	86.55
test1.fcs	20002	1176.59	86.55	86.55
test000.fcs	20004	1176.71	87.73	87.73
test000.fcs	20004	1176.71	87.73	87.73
Gate: Live				
cont1.fcs	18263	2029.22	72.42	72.42
cont1.fcs	18263	2029.22	72.42	91.31
cont1.fcs	18263	2029.22	72.42	72.42
test1.fcs	18328	1078.12	79.30	79.30
test1.fcs	18328	1078.12	79.30	91.63
test1.fcs	18328	1078.12	79.30	79.30
test000.fcs	18868	992.24	73.98	73.98
test000.fcs	18868	992.24	73.98	84.32
test000.fcs	18868	992.24	73.98	73.98

Statistics export

- CSV spreadsheet

Statistics Table




Experiment: **6 color immuno**
Group: **Default_Group_Name**
Sample: **NOT-IN_STATS**

Name	Count	%Gated	%Total
All Events	30000	100.00	100.00
Lymphocytes	6728	22.43	22.43
CD4+	2893	43.00	9.64
CD8+	1358	20.18	4.53
Monocytes	2080	6.93	6.93
Granulocytes	14933	49.78	49.78

- To display **Workspace Statistics Table**, click Statistics without selecting a plot. Workspace statistics contains data of all the gates in the Workspace.
- To display **Plot Statistics Table**, select a plot in the Workspace and then click Statistics. Local statistics only displays data pertaining to the selected plot.
- Alternatively, you can insert Statistics table by right-clicking on a plot or on the workspace, and select Insert Statistics
- Prior to adding a statistics box, make sure the workspace has at least one plot.


Customize Statistics



Experiment: 6 color immuno
Group: Default_Group_Name
Sample: NOT-IN_STATS

Name	Count	%Gated	%Total
All Events	30000	100.00	100.00
Lymphocytes	6728	22.43	22.43
CD4+	2893	43.00	9.64
CD8+	1358	20.18	4.53
Monocytes	2080	6.93	6.93
Granulocytes	14933	49.78	49.78

- To customize Statistics table, select the Table and check statistics to display in the *Statistics* Tab




- To customize Statistics value displayed on a plot, select the plot and choose the statistic

Statistical Values


- **Count:** Number of events collected
- **Events/ μ l:** Concentration of events/ μ l in the gated region
- **% Total:** Percentage of total events collected
- **% Parent:** Percentage of a population based on the number of events collected in the parent gate
- **Mean:** Sum of the signal intensities of a gate divided by the number of values
- **Median** (50th percentile): signal intensity of a gate separating the higher half of a data population
- **Mode:** signal intensity that appears most often in a set of data
- **SD:** Standard Deviation, amount of dispersion of signal intensity around the Mean
- **rSD:** Robust Standard Deviation, amount of dispersion of signal intensity around the Median
- **%CV:** Percent coefficient of variation, Standard Deviation of the peak divided by the Mean of the peak, times 100
- **%rCV:** Percent Robust coefficient of variation, Standard Deviation of the peak divided by the Median of the peak, times 100

Sample Concentration

The **Concentration Statistic** can be selected from the Statistics Ribbon




- Values are displayed as Events/μL



Parameters: BL1-H
Gate: All Events
Experiment: 8 peak beads
Group: Group
Time Recorded: 10:45:51

	Name	Gate	Count	%Total	%Gated	X Median	Concentration
▴	All Events	All Events	10,000	100.000	100.000	26,705	161.3
▴	Beads	Beads	9,038	90.380	90.380	27,084	145.8
▴	R1	R1	1,166	11.660	12.901	467,501	18.8

Customize Panel



Customize

General


Plot Type: ☐ Histogram ☐ Dot ☒ Density ☐ Precedence Density

Resolution: 256 x 256

Color Mode: Log

Color: [Color bar]

% of Events: 100%



Text

Plot: Plot Title 10


Include ☐ Experiment ☐ Group ☒ Sample ☒ Parent Gate

X Axis: FSC-H 9

Y Axis Title: SSC-H

Include ☐ Channel Name ☒ Marker Name ☒ Label Name

Font: T Tahoma



Customize


General

X axis

Y axis

Text

X and Y axis



X axis

Parameter: FSC-H

Scale: ☒ Linear ☐ Logarithmic ☐ HyperLog™

Range: ☐ Automatic ☒ Manual Minimum: 1 Maximum: 1048576

Y axis

Parameter: SSC-H

Scale: ☒ Linear ☐ Logarithmic ☐ HyperLog™


Range: ☐ Automatic ☒ Manual Minimum: 1 Maximum: 1048576

Scale options
Range adjustments

Note: Use < 50 characters when naming gates.
Quad gate names > 31 characters will be truncated.

Heat Map View





- The *Heat Map* view provides a graphical method (i.e., *Heat Map*) for setting up and analyzing plate- and tube-based experiments
 - Access from the “Heat Map” tab of the workspace desktop
 - Displays plate and tube samples



For tubes and plates

Heat Map Setup



- Each sample with a saved FCS file will be colored to reflect a user-specified statistic for a specific gate and parameter
- Select statistic, gate, and parameter for heat map display from drop down menus on the Heat Map Setup Menu



**only gate specific parameters are available*

Heat Map Setup

- Statistical data is overlaid on heat map by selection of the “Show statistic” box under Display mode





- A legend is displayed on the bottom of the Heat Map Display
 - Indicates experiment name, color scheme and transition points

Heat Map Setup

Two display modes are available: Threshold and Heat Map Mode

1. Threshold mode:



- Data displayed using discrete colors to indicate user-specified transition points in data set
- Once level is exceeded, color will change
- Color scheme may be changed by selection of panel of colors



Heat Map Setup

2. Heat Map Mode:



- Data displayed using color gradient to indicate user-specified transition points in data set
- Multiple color choices available from drop down menu



Heat Map Setup

In both modes:

- Define min/max range and transition values by typing value in text boxes
- Add transition points by clicking on colored bar
- Transition points can be repositioned by selecting and dragging arrow to new position
- Delete transition points by clicking on arrow and dragging it away from the colored bar



Set up user-specific filter configurations




Add new filters to the configuration

Add new labels to the detectors/channels

Filter Configuration Manager – Instrument Manager

Allows creation of experiments while not connected to an Attune NxT instrument by maintaining correct instrument settings and channel mapping



Configurations are user-specific – but configurations can be exported and shared .aic files



TIP: the name of the configuration shown is located under the configuration tab

Filter Configuration Manager

The filter labels shown in the configuration display can be viewed by “Label” or by “Wavelength”




Filter Configuration



530/30 as defined in current configuration


Filter Configuration





FITC as defined in current configuration

ADD a new filter to the drop down label menu

1. Instrument > Configuration > Edit Filters
2. Input information about the new filter:
 - Laser Line
 - Filter Type:
 - Collection: bandpass (BP)
 - Directing: Longpass (LP) or Dichroic LP filter
 - Wavelength – specific wavelengths collected or directed by filter



Create new filter or add *New label* to existing filter





Account and Data Management

Attune NxT v 2.5 SW Account Permissions

Permission	Description	User	Advanced User	Administrator	System Administrator
Run Performance Test	Allows a user to run the Performance Test and view the Performance Test reports	X	X	X	X
Run Baseline calculations	Allows a user to run and set new Performance Test Baseline calculations		X	X	
Advanced instrument settings	Allows a user to adjust the width threshold, window extensions, and area scaling factor within an experiment		X	X	X
Run system decontamination	Allows a user to run the Decontaminate System function		X	X	
Run System Tests	Allows a user to run the system tests (Button on instrument ribbon is labeled "Self Test")		X	X	
Manage User accounts	Allows a user to create user accounts, edit user accounts, reset passwords, change passwords, view login/logout times for all users, and view the length of all user sessions			X	X
Set security policy	Allows a user to set system security settings for username length, password length, password expiration and lock-out, and the auto lock out time due to system inactivity				X



System Logs / User Logs



- **General Log**
- Log of system transactions.
- The information displayed for administrator
- Important for troubleshooting/support
- **User Log**
- Filter by:
 - Date
 - User
 - Elapsed time or Sample Count
- Shared instrument? – helpful for tracking usage

Export Statistics to CSV


- Select the Experiment, Group, or Sample files to export
- Multi-selection enabled with the Ctrl button
- Right-click and select Export
- Select statistics level and views
- Select single or separate files



Results table

- View the statistics associated in a tabulated format


Select statistics



Group results


Filename	Gate	Count	Concentration	%Total	%Gated
Gate: All Events					
wt.fcs	All Events	218...	2182.54	100.00	100.00
wt.fcs	All Events	218...	2182.54	100.00	100.00
wt.fcs	All Events	218...	2182.54	100.00	100.00
wt.fcs	All Events	218...	2182.54	100.00	100.00
KI.fcs	All Events	141...	2184.29	100.00	100.00
KI.fcs	All Events	141...	2184.29	100.00	100.00
KI.fcs	All Events	141...	2184.29	100.00	100.00
KI 200.fcs	All Events	144...	2217.83	100.00	100.00
KI 200.fcs	All Events	144...	2217.83	100.00	100.00
KI 200.fcs	All Events	144...	2217.83	100.00	100.00
Gate: R1					
wt.fcs	R1	174...	1749.83	80.17	80.17
wt.fcs	R1	174...	1749.83	80.17	80.17
KI.fcs	R1	1135...	1747.57	80.01	80.01
KI.fcs	R1	1135...	1747.57	80.01	80.01
KI 200.fcs	R1	1140...	1755.29	79.14	79.14
KI 200.fcs	R1	1140...	1755.29	79.14	79.14
Gate: R2					
wt.fcs	R2	3373	33.73	1.55	1.93
wt.fcs	R2	3373	33.73	1.55	1.55
wt.fcs	R2	3373	33.73	1.55	1.55
KI.fcs	R2	5079	78.14	3.58	4.47
KI.fcs	R2	5079	78.14	3.58	3.58
KI.fcs	R2	5079	78.14	3.58	3.58
KI 200.fcs	R2	5096	78.40	3.53	4.47
KI 200.fcs	R2	5096	78.40	3.53	3.53
KI 200.fcs	R2	5096	78.40	3.53	3.53

Export



Print Results



- Select the Experiment, Group, or Sample files
- Multi-selection enabled with the Ctrl or shift button
- Select workspace level and views
- Select print destination



Experiment Export


- Right click on the experiment name, open context menu
- Select *Export Experiment*
- Windows browser opens to last 'saved' directory
- Creates files:

.atx tube experiment
.apx plate experiment




Export/Save of files as FCS 3.0 and FCS 3.1


All experiment samples



Selected experiment samples



Right click on the experiment name or selected sample(s) and select “Export FCS file” from the drop down menu




FCS 3.1 files are compatible with: FlowJo V10+, Kaluza, FCS Express V4+, and VenturiOne 3rd

FCS 3.0 files are compatible with: FlowJo V7+, FCS Express V3+, VenturiOne, and Kaluza

Experiment Import


Right click on the *User* name, open context menu select *Import Experiment*



Windows browser opens to last 'saved' directory

Select:

.atx tube experiment
.apx plate experiment



File Types and Extensions

Data format: FCS 3.1 (default)

Storage location: the directory where last file was saved

File extensions - automatically added to each file

.fcs – data file

.ahm – heat map file

.arp – run protocol*

.aws – workspace file*

.ais – instrument settings file*

.acs – compensation settings*

.aic – instrument configuration*

.afs – system log in system local format



* These files can be saved (exported) and used (loaded) in other experiments

If samples are named using any of the following words, the FCS file will not be recorded:



CON, PRN, AUX, CLOCK\$, NUL, COM1, COM2, COM3, COM4, COM5, COM6, COM7, COM8, COM9, LPT1, LPT2, LPT3, LPT4, LPT5, LPT6, LPT7, LPT8, LPT9

Templates

- Experiment Templates are directions for a new experiment that include workspace objects, run protocol, and instrument settings as defined by the user
- Templates can be accessed from a single user profile
- To create a template from a recorded experiment, select “Save as Template” from Home Tab or by right-clicking on user name in Experiment Explorer
- To create a new experiment from template, select the “Templates” button located on the Main Menu desktop or from the Home Tab
- Once a Template is created, it cannot be deleted



System Log



User Log

Filter by:


- Date
- User
- Elapsed time or Sample Count

General Log

- Log of system transactions.
- The information displayed for administrator

Open the *View* tab

Check the FCS file panel



FCS File

- Floating panel with collapsible sections
 - File name and path
- FCS file info
- Sample information
 - Start/end time
 - Flow rate
 - Volume
 - Total events
 - Lost events
- Parameters
 - Channel
 - Target & label
 - voltage
- Compensation
 - Spillover values
- System information
 - Configuration
 - Laser, ASF, laser delay

The screenshot shows a software window titled "FCS Information" with a dark theme. It contains three main sections: "FCS file information", "Sample Information", and "Parameters (13)".

FCS file information

Filename: 4 COLOR.fcs
Filepath: C:\Users\Public\Documents\Life Technologies\Attune\

Sample Information

Filename (\$FIL) Sample.fcs
Date (\$DATE) 17-Mar-2015
Start Time (\$BTIM) 15:31:12
End Time (\$ETIM) 15:31:42
Institution (\$INST) NA
Investigator (\$EXP) NA
operator (\$OP) NA
Originality (\$ORIGINA... Original
Experiment (\$PROJ) concentration
Specimen (\$SMNO) LYSED WHOLE BLOOD
Sample Source (\$SRC) NA
Sample Information (\$... NA
Flow Rate (#FLOWRA... 100
Comment (\$COM) NA
Volume (\$VOL) 50000
Timestep (\$TIMESTEP) 0.001
Trigger 1 (#TR1) FSC,150000
Lost Events (\$LOST) 0
Total Events (\$TOT) 16191

Parameters (13)

Param#	1	2	3	4	5
Name	Time	FSC-A	SSC-A	BL3-A	YL
Stain	Time	FSC-A	SSC-A	cd19-PER...	CD3
Target	NA	-	-	cd19	CD3
Label	NA	-	-	PERCP-C...	PE
Voltage	NA	370	400	450	310
Range	67108864	1048576	1048576	1048576	104
Bits	32	32	32	32	32
Excitation Wavelength...	NA	NA	NA	NA	NA
Optical Filter	NA	NA	NA	NA	NA
Amplification	0,0	0,0	0,0	0,0	0,0

© 2020 Thermo Fisher Scientific, Inc. All rights reserved.

All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For Research Use Only. Not for use in diagnostic procedures.